

# Journal Pre-proof

"How I do it": Home Non-Invasive Ventilation in COPD

Marta Kaminska, MD, Veronique Adam, RRT, Jeremy Orr, MD

PII: S0012-3692(24)00041-2

DOI: <https://doi.org/10.1016/j.chest.2024.01.030>

Reference: CHEST 6059

To appear in: *CHEST*

Received Date: 25 September 2023

Revised Date: 19 December 2023

Accepted Date: 15 January 2024

Please cite this article as: Kaminska M, Adam V, Orr J, "How I do it": Home Non-Invasive Ventilation in COPD, *CHEST* (2024), doi: <https://doi.org/10.1016/j.chest.2024.01.030>.

This is a PDF file of an article that has undergone enhancements after acceptance, such as the addition of a cover page and metadata, and formatting for readability, but it is not yet the definitive version of record. This version will undergo additional copyediting, typesetting and review before it is published in its final form, but we are providing this version to give early visibility of the article. Please note that, during the production process, errors may be discovered which could affect the content, and all legal disclaimers that apply to the journal pertain.

Copyright © 2024 Published by Elsevier Inc under license from the American College of Chest Physicians.



**Title:** "How I do it": Home Non-Invasive Ventilation in COPD

**Authors:** Marta Kaminska<sup>1,2</sup> , MD, Veronique Adam RRT<sup>1</sup> , Jeremy Orr, MD<sup>3</sup>

<sup>1</sup>Quebec National Program for Home Ventilatory Assistance-McGill University Health Center (PNAVD-MUHC), Montreal, QC, Canada

<sup>2</sup>Division of Respiratory Medicine, Department of Medicine, McGill University Health Centre, Montreal, QC, Canada

<sup>3</sup>University of California, Division of Pulmonary, Critical Care, Sleep Medicine, and Physiology, La Jolla, California, USA.

**Corresponding Author:**

Marta Kaminska, Respiratory Epidemiology and Clinical Research Unit, Centre for Outcomes Research and Evaluation, Research Institute of the McGill University Health Centre, 5252 de Maisonneuve Boulevard West, Montreal, QC, H4A 3S9, Canada. Email: [marta.kaminska@mcgill.ca](mailto:marta.kaminska@mcgill.ca)

**Key Words:**

Chronic Obstructive Pulmonary Disease; Chronic Respiratory Failure; Non-Invasive Ventilation; Hypoventilation

**Disclosures:**

Dr. Kaminska reports research grants from the CHEST foundation and CIHR, investigator-initiated research support from Philips, and advisory board membership for Biron Soins du Sommeil. Veronique Adam: none. Dr. Orr participates in an advisory board for Resmed Inc. and receives research funding from NIH/NHLBI.

**Acknowledgement:**

We would like to thank the staff of the Quebec National Program for Home Ventilatory Assistance-McGill University Health Center ("PNAVD") for their dedication to their patients and the excellence of their work.

**Word Count:**

ABSTRACT 163

BODY 3193

**Abbreviations**

ATS	American Thoracic Society
AVAPS/iVAPS	Average / Intelligent Volume Assured Pressure Support
BMI	Body Mass Index
BUR	Back-up Rate
COPD	Chronic Obstructive Pulmonary Disease
CTS	Canadian Thoracic Society
EPAP	Expiratory Positive Airway Pressure
ERS	European Respiratory Society
EFL	Expiratory flow limitation
FEV <sub>1</sub>	Forced expiratory volume in 1 second
IPAP	Inspiratory Positive Airway Pressure
LT-NIV	Long-term Non-Invasive Ventilation
LTOT	Long-term Oxygen Therapy
MPV	Mouthpiece Ventilation
NIV	Non-Invasive Ventilation
PCO <sub>2</sub>	Partial pressure of carbon dioxide
PEEP	Positive end-expiratory pressure
PS	Pressure Support
RT	Respiratory Therapist
TcCO <sub>2</sub>	Transcutaneous carbon dioxide
Ti	Inspiratory time
QoL	Quality of life

**Abstract**

There is increasing evidence that long-term NIV (LT-NIV) can improve outcomes in individuals with severe, hypercapnic COPD. Though the evidence remains unclear in some aspects, LT-NIV appears to be able to improve patient-related and physiological outcomes like dyspnea, FEV<sub>1</sub>, PCO<sub>2</sub>, and also reduce rehospitalizations and mortality. Efficacy is generally associated with reduction in PCO<sub>2</sub>. To achieve this, an adequate interface (mask) is essential, as are appropriate ventilation settings that target the specific respiratory physiology of COPD. This will ensure comfort, synchrony and adherence that will result in physiologic improvements. This article briefly reviews the newest evidence and current guidelines on LT-NIV in severe COPD. It describes a true case who benefitted from the therapy. Finally, it provides strategies for initiating and optimizing this LT-NIV in COPD, discussing high-pressure NIV, optimization of triggering and control of inspiratory time. As demand increases, clinicians will need to be familiar with this therapy, to reap its benefits, as inadequately adjusted LT-NIV will not be tolerated or effective.

## 1. Introduction

COPD is a major cause of morbidity and the third-leading cause of death worldwide<sup>1</sup>. With progressive disease, chronic hypercapnic respiratory failure may occur, which is associated with a 1-year mortality of 17-30%<sup>2-5</sup>. Hospitalization due to hypercapnic exacerbation is associated with 11% in-hospital mortality and 49% 2-year mortality<sup>6</sup>.

The benefits of non-invasive ventilation (NIV) in acute exacerbations are well-established<sup>7-10</sup>. For chronic use in the home setting, older systematic reviews found no evidence of benefit but identified higher baseline CO<sub>2</sub>, higher pressures and longer daily NIV use as predictors of greater CO<sub>2</sub> reduction<sup>11</sup>. Recent randomized controlled trials have changed the landscape in the field. However, practice had been variable, from higher proportions of patients with COPD in ventilation programs in European countries<sup>12-14</sup> to few across Canadian provinces<sup>15</sup>, and generally few chronically hypercapnic COPD patients receiving home NIV in the United States<sup>16</sup>. Hence, experience with long-term NIV (LT-NIV) is highly heterogeneous among pulmonologists. Varying ventilation targets<sup>17</sup>, rapid technological advances with different ventilation modes, diverse settings and nomenclature between manufacturers, and the multitude of parameters on NIV devices resulting in wide range of possible prescriptions, all make NIV a complex therapy that may be challenging to implement.

The purpose of this review is to summarize the current data and guidelines on LT-NIV in COPD and provide a practical approach to its initiation.

## 2. Evidence for LT-NIV in COPD

Randomized controlled trials of NIV in COPD published prior to 2010 largely failed to demonstrate benefits on physiologic, functional or patient-reported outcomes, reduction of hospitalizations or survival<sup>11,18</sup>. One exception was an Australian trial comparing LT-NIV with long-term oxygen therapy (LTOT) or LTOT alone which found improvement in mortality with LT-NIV after a mean follow-up of 2.21 years<sup>4</sup>. However, there was worsening of certain quality of life (QoL) parameters, which hampered widespread support for LT-NIV in COPD<sup>19</sup>.

In parallel, work done in Europe generated increasing evidence that LT-NIV could result in improved outcomes, if settings were adjusted targeting obstructive COPD physiology and with specific CO<sub>2</sub> reduction goals, termed “high-intensity” ventilation<sup>20-22</sup>. While the exact definition of high-intensity NIV remains somewhat contentious, it is primarily the high pressure rather than

a high back-up rate that is the most relevant parameter<sup>23,24</sup>. Using this type of strategy, Kohnlein et al. demonstrated markedly improved survival with LT-NIV in stable hypercapnic COPD patients<sup>5</sup>, with improved QoL. The mean inspiratory (IPAP) and expiratory (EPAP) pressures were 21.6 and 4.8 cmH<sub>2</sub>O respectively, with back-up rate (BUR) 16.1/min. This trial included participants with PaCO<sub>2</sub>  $\geq$  52 mmHg (mean 58 mmHg) and showed reduction, but not normalization of daytime PaCO<sub>2</sub>, and improvement in oxygen saturation and FEV<sub>1</sub>. Murphy et al. studied patients who remained hypercapnic (mean PCO<sub>2</sub> 59 mmHg) 2-4 weeks following hospitalization for exacerbation requiring acute NIV, and who met criteria for LTOT<sup>25</sup>. Using mean IPAP and EPAP of 24 and 4 cmH<sub>2</sub>O, respectively, they demonstrated delayed time to readmission, reduced readmission rate and improved QoL. There was no improvement in mortality but many control patients eventually crossed over to active treatment, and there were few deaths. Conversely, another trial also including hypercapnic (mean PCO<sub>2</sub> approximately 59 mmHg) participants post exacerbation did not demonstrate benefit of NIV on readmissions or death<sup>26</sup>, using mean IPAP and EPAP 19 and 4 cmH<sub>2</sub>O, respectively. It has been suggested this resulted from inclusion of participants randomized as early as 48 hours after discontinuation of acute NIV, who may not have been chronically hypercapnic, since reduction of PCO<sub>2</sub> in the control group mirrored that in the NIV group in the first 3 months. Variable results across modern trials highlight that LT-NIV benefits very specific COPD patient populations; the positive trials define the indications and target populations.

Real world data, notwithstanding potential biases, have helped understand the impact of LT-NIV in North America. Retrospective single center studies demonstrated good adherence<sup>27</sup>, improved event-free survival with LT-NIV post discharge<sup>28</sup> and marked reduction in re-hospitalizations with a multimodal intervention including LT-NIV<sup>29</sup>. Using Medicare data, it was demonstrated that in patients with COPD and chronic respiratory failure, LT-NIV was associated with reduced all-cause mortality, hospitalizations, and ER visits<sup>30,31</sup>, with benefit on mortality and reduction in Medicare expenditures restricted to those with hypercapnia<sup>16</sup>.

### 3. Clinical Practice Guidelines

The European Respiratory Society (ERS)<sup>32</sup>, American Thoracic Society (ATS)<sup>33</sup> and the Canadian Thoracic Society (CTS)<sup>34</sup> have recently published guidelines that agree there is a role for LT-NIV in chronic stable hypercapnic COPD, but recommendations are weak, with very low

to moderate certainty evidence. LT-NIV is suggested following acute hypercapnic respiratory failure requiring acute NIV, provided the patient remains hypercapnic after resolution of the acute episode. The ATS and CTS recommend re-evaluation after 2-4 weeks. All three guidelines specify that NIV should aim to normalize or at least significantly reduce hypercapnia. The ERS and CTS found no evidence to preferentially support auto-adjusting modes and recommend fixed pressure modes. The ATS suggest evaluating patients for obstructive sleep apnea (OSA), at a minimum with a questionnaire (very low certainty evidence). There remain many unanswered questions regarding implementation of LT-NIV in COPD, including the need to better identify patients most likely to benefit, develop strategies to initiate and optimize LT-NIV, and clarify management of obese hypercapnic COPD patients, especially given high obesity rates in North America including in COPD<sup>27</sup>.

#### 4. Case Study

A 57-year-old male was referred for optimization of home NIV in December 2021. He was a former smoker of 35 pack-years, who quit 1 year prior, with body mass index (BMI) 25 kg/m<sup>2</sup>, FEV<sub>1</sub> 0.6L (19% predicted), on optimal bronchodilator therapy and LTOT. He had multiple hospitalizations for hypercapnic respiratory failure over the 2 years prior, with maximal venous PvCO<sub>2</sub> 136 mmHg, requiring NIV acutely several times and intubation once. His arterial blood gas 3 years prior showed pH 7.4 and PaCO<sub>2</sub> 59 mmHg.

*This patient would have likely benefitted from earlier implementation of home NIV, which may have prevented re-hospitalizations and intubation.*

He had been started on nocturnal NIV at the referring center using an oronasal mask, in spontaneous-timed (ST) mode as per parameters used during his last admission: IPAP 16 cmH<sub>2</sub>O, EPAP 10 cmH<sub>2</sub>O, BUR 12/min, rise time 300 ms, inspiratory time (Ti) minimum 0.8 sec and maximum 1.6 sec, high trigger sensitivity and medium cycling sensitivity.

*These parameters are not optimal due to low IPAP and low driving pressure ( $\Delta$ =IPAP-EPAP). It is unclear why the EPAP was relatively high but this was comfortable for the patient. The rise time should be faster, Ti maximum lower and cycling sensitivity high for earlier cycling from IPAP to EPAP (cf below).*

Despite symptomatic improvement, he was briefly re-hospitalized with  $P_vCO_2$  102 mmHg. Upon discharge, parameters were modified remotely by his home care respiratory therapist (RT) to IPAP 19 cmH<sub>2</sub>O, with fastest rise time. The patient reported good subjective sleep quality with NIV. A few days later, a routine follow-up  $P_vCO_2$  taken at home was 119 mmHg.

*The driving pressure remained low in the context of COPD. Note that  $P_vCO_2$ , while convenient, is not recommended as it is not a reliable estimate of  $P_aCO_2$ .*

IPAP was gradually increased remotely to 25 cmH<sub>2</sub>O, with Ti max 1.3 sec. About 9 months after starting NIV, an overnight oximetry on NIV and O<sub>2</sub> at 3L/ minute showed mean SpO<sub>2</sub> 95% ,and SpO<sub>2</sub> <90% during 0.8% of the night. His daytime transcutaneous CO<sub>2</sub> (TcCO<sub>2</sub>) taken in the outpatient clinic was 58 mmHg. Supplemental oxygen at night was reduced.

*Excessive O<sub>2</sub> supplementation may result in CO<sub>2</sub> retention. O<sub>2</sub> should be titrated to maintain SpO<sub>2</sub> not higher than 92% at rest.*

EPAP was decreased to 8cmH<sub>2</sub>O and Ti range to 0.8 – 1.1 sec, but an attempt at increasing IPAP was not tolerated. Overnight TcCO<sub>2</sub> showed a mean of 54mmHg, fluctuating depending on leak (incompletely controlled). Due to increasing daytime NIV use, mouthpiece ventilation (MPV) was added to help with dyspnea during daily activities. It allowed the patient to do some work in his garage, for which he was grateful. About 24 months after starting NIV, he remains at home with no further hospitalizations despite exacerbations requiring antibiotic and corticosteroid treatment.

*Despite imperfect CO<sub>2</sub> control, the patient was stabilized and re-hospitalizations have to date been averted.*

## 5. Approach to home NIV in COPD

### 5.1. Initiation setting

There is no agreed-upon optimal initiation location. The ATS suggests *against* titrating NIV in the sleep laboratory due to concerns regarding cost, delay, safety of achieving rapid normocapnia over a single night and proficiency of personnel in such titrations<sup>33</sup>. European studies used initiation in hospital, often in specialized ventilation units, over days, with progressive acclimatization and optimization<sup>5,17,25</sup>. This is not currently possible in most North



American centers. However, parameters used and tolerated during an acute exacerbation can be a useful starting point. Titration performed while admitted for an acute exacerbation is practical, but chronic hypercapnia should be confirmed before starting LT-NIV and parameters reassessed based on ventilation goals.

Home initiation has been studied as an alternative to in-hospital initiation<sup>35</sup>. The protocol required multiple home visits and daily calls from a specialized nurse, remote monitoring of ventilator and overnight TcCO<sub>2</sub> data, and lasted a median of 7 days. Technical difficulties affected TcCO<sub>2</sub> measurements in a significant proportion of patients. This protocol was non-inferior to hospital initiation, cheaper and safe (patients with unstable cardiac comorbidities or heart failure were excluded). This highly resource-intensive protocol is promising but would likely require further optimization and adaptation to local realities.

In our experience, NIV initiation often takes place in the outpatient setting, utilizing a combination of resources, ideally specialized RTs, available in clinic and at home. Clinic visits may include discussions regarding NIV potential benefits and challenges with shared decision-making regarding initiation, training on equipment, and interface selection. A daytime NIV trial may be performed, starting at low pressures and increasing pressures as tolerated during quiet wakefulness. Alternatively, initiation may be performed directly at home, utilizing trained home care or medical equipment providers.

Close follow-up in the initial period after starting NIV is important. Trained staff (ideally an RT) should follow-up within days to provide support and optimize interface and ventilation comfort. This could be done remotely, including for parameter adjustments (for non-life support devices)<sup>35</sup>, providing there is adequate connectivity. New mask trials, if needed, should be done in person. The extent to which this is done by home care or equipment companies versus clinic staff may vary based on local resources. Physician follow-up should occur within weeks to assess adherence, comfort and PCO<sub>2</sub> [arterial, capillary or TcCO<sub>2</sub>] to inform further parameter adjustments. We suggest targeting reduction of daytime PCO<sub>2</sub> by  $\geq 20\%$ , or to  $< 48$  mmHg if possible<sup>5</sup>. Adjustments can then be completed targeting nocturnal gas exchange<sup>25</sup> based on overnight TcCO<sub>2</sub> or oximetry, remembering the latter will not be reliable regarding hypoventilation in patients on LTOT. This type of outpatient protocol may take longer to reach optimal parameters than hospital-based initiation but may constitute a lesser burden on patients and healthcare resources.

## 5.2. Mask choice:

The choice of mask is crucial for initiation and adaptation to LT-NIV, and is a key issue for the majority of NIV users and caregivers<sup>36</sup>. Oro-nasal interfaces appear to be most frequently used<sup>36</sup>, possibly because of oro-nasal mask use during hospitalizations. Additionally, very dyspneic patients may be unable to breathe exclusively through their nose, precluding use of a nasal interface. However, nasal masks should still be considered, as they cover a smaller facial area, hence are less prone to leaks, and might be easier to apply. A chin strap can help counter mouth leak<sup>37</sup> for patients able to use one. Nasal masks also have less potential for exacerbating upper airway obstruction<sup>38,39</sup> but this has not been studied in the context of NIV.

To improve adherence, the chosen mask must be comfortable and easy to handle. Patients applying their mask independently must be able to do it easily and quickly if dyspneic. If experiencing difficulty lifting arms up to the head, strategies should be developed to limit required movement during mask application. A well-fitting mask can determine adherence by minimizing uncomfortable pressure points, and preventing leaks which can be bothersome, cause asynchrony, reduce effectiveness of ventilation and introduce errors into reported ventilation.

## 5.3. NIV Parameters

### Pressures

Achieving significant PCO<sub>2</sub> reduction is the main goal of NIV, along with resting of the diaphragm and accessory muscles. Using bilevel NIV, this typically requires high driving pressures (IPAP >18 cm H<sub>2</sub>O) to overcome high airways resistance present also in inspiration in severe COPD. This results in larger tidal volumes, allowing adequate alveolar ventilation in the context of increased dead space, and lower respiratory rates to help reduce hyperinflation<sup>22</sup>. In a naïve patient, IPAP can be started at 12 cmH<sub>2</sub>O, and increased empirically to maximal tolerated, aiming a normalizing PCO<sub>2</sub><sup>22</sup>. Pressures up to 42 cmH<sub>2</sub>O have been reported<sup>40</sup>, though values above 30 cmH<sub>2</sub>O are unusual. While IPAP may need to be gradually increased to allow patient adaptation, it has been shown that adherence is better with higher compared to lower pressures, in addition to improved PCO<sub>2</sub>, dyspnea and FEV<sub>1</sub><sup>21</sup>.

Volume assured pressure modes (AVAPS/iVAPS), with or without auto-EPAP, can be a good option for auto-titration or for long-term use. However, mask leaks must be controlled for device algorithms to function properly. They appear to be equivalent with respect to efficacy to fixed PS modes<sup>41,42</sup>, though some evidence suggests auto-EPAP can be superior regarding PCO<sub>2</sub> control and symptoms<sup>43,44</sup>. A starting target tidal volume of 6-8 mL per kilogram of ideal body weight can be used, with subsequent adjustment to optimize PCO<sub>2</sub> and comfort. The set minimal IPAP should be sufficiently high to ensure the target volume is reached through support provided primarily by the device rather than patient-generated effort.

The EPAP has traditionally been set at low values in non-obese COPD patients. A level of 3-5 cmH<sub>2</sub>O should be used to start (devices limit minimum EPAP of 3-4 cmH<sub>2</sub>O to prevent re-breathing in a single limb passive circuit). However, higher EPAP may be needed to correct upper airway obstruction<sup>45</sup>, and can be titrated during polysomnography or using an auto-EPAP mode. EPAP can also help match intrinsic positive end-expiratory pressure (auto-PEEP), reducing triggering effort and asynchrony<sup>46</sup>. However, auto-PEEP cannot easily be measured. There is some evidence that higher EPAP may also correct expiratory flow limitation (EFL) by “stenting” the airways during expiration<sup>46</sup> which could improve comfort and reduce hyperinflation. NIV devices with self-adjusting EPAP targeting EFL exist, but are not widely available, and require further study. Manual adjustment optimizing patient device triggering and comfort remains necessary.

#### Back-up rate & Trigger sensitivity

Effective triggering is important for optimal synchrony, to minimize work of breathing and for comfort<sup>12</sup>. In addition to optimizing mask fit and EPAP, this can be achieved by adjusting trigger sensitivity in many devices, starting with medium. Initial high-intensity protocols included high BUR to control ventilation and minimize patient effort<sup>47</sup>. However, it has been shown that with adequate pressures, a high BUR does not provide further benefit or unloading of respiratory muscles<sup>23</sup>. In the US, obtaining a device for COPD with a BUR often requires either prescribing a ventilator, or “failing” a spontaneous bilevel device. For less severe patients, starting without a BUR may be reasonable. For more severe patients, obtaining a device with a BUR is advisable. When utilizing a BUR, we suggest starting at 12/min and titrating up only after optimizing

pressures, triggering, and Ti/cycling. Overly aggressive BUR may exacerbate patient-ventilator asynchrony (e.g. glottic closure)<sup>48</sup> and promote intolerance.

#### Ti and Cycling Off

Longer expiratory time helps reduce hyperinflation in COPD, hence a shorter Ti (Figure 1) is beneficial provided ventilation remains adequate. Delayed cycling off (returning from IPAP to EPAP)<sup>49</sup> can lead to asynchrony during flow-based cycling modes like bilevel S or ST, and VAPS (Figure 2). Cycling off should therefore be set at high sensitivity (earlier). Uncontrolled leaks may prevent synchronized cycling. Pressure (assist) control (PC / PAC) modes, where the Ti is fixed for all breaths, can be useful to remedy these problems. For bilevel S or ST and some VAPS modes, many devices allow limits on Ti on spontaneous breaths (min/max Ti); device-triggered breaths have a fixed Ti. We suggest a Ti min and max starting range of 0.5 - 1.0 sec.

#### Rise time

Patients with COPD often have “air hunger”. A rapid pressure rise from EPAP to IPAP provides high flow when inspiration is triggered, reducing work of breathing<sup>50</sup>. This also allows reaching target IPAP faster and higher delivered volume for a given Ti.

#### Ramp up and Ramp down:

Bilevel devices may have a “Ramp up” feature allowing pressures to increase progressively when started. This can help when high IPAP is difficult to tolerate. Yet if pressures are too low, this may be uncomfortable and result in uncontrolled hypoventilation<sup>51</sup>, hence this should be applied judiciously and followed-up.

Some COPD patients experience prolonged dyspnea upon stopping NIV in the morning, termed deventilation syndrome<sup>52</sup>. It can be remedied by adjusting parameters and improving synchrony<sup>52</sup>. Some devices now also include a “Ramp down” feature where pressures decrease gradually for an intended smoother transition to spontaneous breathing, but no clinical data are available on this feature.

#### 5.4.Comorbidities

A significant proportion of COPD patients also have OSA, obesity, or another restrictive condition such as kyphoscoliosis. OSA can be corrected by adapting EPAP. However, the literature is scant on how to approach other parameters in the context of mixed obstructive and restrictive syndromes. While NIV for restrictive conditions also requires high IPAP, one would typically use longer T<sub>I</sub>, slower rise time and delayed cycling off. Our practice is typically dictated by the primary spirometric abnormality. In a patient with milder COPD with a primarily restrictive syndrome, we chose parameters targeting the restriction. Adjustments are made as needed to optimize comfort, synchrony and gas exchange.

### 5.5. Monitoring

Monitoring of symptoms and daytime PCO<sub>2</sub> are key elements providing the clinician with feedback regarding NIV effectiveness. Overnight gas exchange should be assessed<sup>53</sup>. Oximetry is frequently used but difficult to interpret with supplemental oxygen. Overnight home TcCO<sub>2</sub> is useful<sup>54</sup> but remains technically challenging, expensive, and not widely available in North America. Device recordings (either via manual “download” or cloud-based remote monitoring) can help ascertain adherence, interface adjustment (leak), and more detailed ventilation outcomes<sup>55</sup>. Importantly, download data, such as tidal volume and minute ventilation, are inaccurate in the face of significant leak.

### 6. Applications beyond night-time use

Very dyspneic patients may benefit from NIV extension into the daytime. This can be achieved using the same settings as for sleep, e.g. during a rest period, but favoring a nasal mask to allow easier speech. Open circuit mouthpiece ventilation (MPV)<sup>56</sup> allows more freedom, mobility and improves ability to communicate. Commercial home ventilators now have dedicated MPV modes and circuits. Hand-held portable devices have been shown to improve dyspnea and walking distance<sup>57</sup>, but can be heavy and impractical.

### 7. Conclusions

Data are accumulating, suggesting home LT-NIV in patients with hypercapnic severe COPD can improve clinical and patient-related outcomes, when the parameters are adapted to COPD physiology. Many questions remain regarding patient selection and implementation. COPD

310 patients bring a different set of challenges compared with other populations requiring home NIV.  
311 Pulmonologists will inevitably be faced with increasing demand that will require locally adapted  
312 protocols and care pathways to optimize delivery of this therapy.

313

Journal Pre-proof

**Table:**

Table. Principles of mask selection and application

Always confirm mask fit with NIV on
Lift the mask to inflate the air cushion, before tightening the headgears
Use the minimum tension required on the headgear
Always try the mask in the position it will be used
Verify that the patient can apply the mask correctly and without excessive effort
Clean the mask daily to remove the film left by skin sebum

**Figure legends:**

Figure 1. Main Elements of a Bi-level NIV Device -Supported Breath.

From Selim BJ, Wolfe L, Coleman JM, 3rd, Dewan NA. Initiation of Noninvasive Ventilation for Sleep Related Hypoventilation Disorders: Advanced Modes and Devices. *Chest* 2018; 153: 251-265.


Figure 2. Cycling threshold affects the duration of respiration. A: Normal respiratory mechanics. The cycle setting is 25% of peak inspiratory flow. Cycling is ideal, as indicated by the fact that the inspiratory flow decreases to the 25% cycling level at the end of the patient's neural inspiratory time ( $T_I$ ). B: Obstructive respiratory mechanics. The change in the inspiratory flow curve leads to the 25% level being reached later, well after the end of the neural  $T_I$ . The duration of delayed cycling is represented by the excess  $T_I$ . Increasing the cycle setting to 60% of peak inspiratory flow corrects this problem, and cycling occurs at the end of neural  $T_I$ . Insp = inspiration. Exp = expiration.

From Gentile MA. Cycling of the mechanical ventilator breath. *Respir Care*. 2011;56(1):52-60

334


335 Figure1 . (Permission will be required: Selim et al. *Chest* 2018; 153: 251)

336





337 Figure 2, (Permission will be required; From Gentile, Respir Care 2011)



338

## References

1. The top 10 causes of death. 2020; <https://www.who.int/news-room/fact-sheets/detail/the-top-10-causes-of-death>. Accessed 11/15/2023, 2023.
2. Casanova C, Celli BR, Tost L, et al. Long-term controlled trial of nocturnal nasal positive pressure ventilation in patients with severe COPD. *Chest*. 2000;118(6):1582-1590.
3. Clini E, Sturani C, Rossi A, et al. The Italian multicentre study on noninvasive ventilation in chronic obstructive pulmonary disease patients. *The European respiratory journal*. 2002;20(3):529-538.
4. McEvoy RD, Pierce RJ, Hillman D, et al. Nocturnal non-invasive nasal ventilation in stable hypercapnic COPD: a randomised controlled trial. *Thorax*. 2009;64(7):561-566.
5. Kohnlein T, Windisch W, Kohler D, et al. Non-invasive positive pressure ventilation for the treatment of severe stable chronic obstructive pulmonary disease: a prospective, multicentre, randomised, controlled clinical trial. *The Lancet. Respiratory medicine*. 2014;2(9):698-705.
6. Connors AF, Jr., Dawson NV, Thomas C, et al. Outcomes following acute exacerbation of severe chronic obstructive lung disease. The SUPPORT investigators (Study to Understand Prognoses and Preferences for Outcomes and Risks of Treatments). *American journal of respiratory and critical care medicine*. 1996;154(4 Pt 1):959-967.
7. Rochweg B, Brochard L, Elliott MW, et al. Official ERS/ATS clinical practice guidelines: noninvasive ventilation for acute respiratory failure. *The European respiratory journal*. 2017;50(2).
8. Chandra D, Stamm JA, Taylor B, et al. Outcomes of noninvasive ventilation for acute exacerbations of chronic obstructive pulmonary disease in the United States, 1998-2008. *American journal of respiratory and critical care medicine*. 2012;185(2):152-159.
9. Elliott MW, Nava S. Noninvasive ventilation for acute exacerbations of chronic obstructive pulmonary disease: "Don't think twice, it's alright!". *American journal of respiratory and critical care medicine*. 2012;185(2):121-123.
10. Lindenauer PK, Stefan MS, Shieh MS, Pekow PS, Rothberg MB, Hill NS. Outcomes associated with invasive and noninvasive ventilation among patients hospitalized with exacerbations of chronic obstructive pulmonary disease. *JAMA internal medicine*. 2014;174(12):1982-1993.
11. Struik FM, Lacasse Y, Goldstein RS, Kerstjens HA, Wijkstra PJ. Nocturnal noninvasive positive pressure ventilation in stable COPD: a systematic review and individual patient data meta-analysis. *Respiratory medicine*. 2014;108(2):329-337.
12. Lloyd-Owen SJ, Donaldson GC, Ambrosino N, et al. Patterns of home mechanical ventilation use in Europe: results from the Eurovent survey. *The European respiratory journal*. 2005;25(6):1025-1031.
13. Cantero C, Adler D, Pasquina P, et al. Long-Term Noninvasive Ventilation in the Geneva Lake Area: Indications, Prevalence, and Modalities. *Chest*. 2020;158(1):279-291.
14. Patout M, Lhuillier E, Kaltsakas G, et al. Long-term survival following initiation of home non-invasive ventilation: a European study. *Thorax*. 2020;75(11):965-973.
15. Rose L, McKim DA, Katz SL, et al. Home mechanical ventilation in Canada: a national survey. *Respir Care*. 2015;60(5):695-704.

16. Frazier WD, DaVanzo JE, Dobson A, Heath S, Mati K. Early Initiation of non-invasive ventilation at home improves survival and reduces healthcare costs in COPD patients with chronic hypercapnic respiratory failure: A retrospective cohort study. *Respiratory medicine*. 2022;200:106920.
17. van der Leest S, Duiverman ML. High-intensity non-invasive ventilation in stable hypercapnic COPD: Evidence of efficacy and practical advice. *Respirology (Carlton, Vic.)*. 2019;24(4):318-328.
18. Struik FM, Lacasse Y, Goldstein R, Kerstjens HM, Wijkstra PJ. Nocturnal non-invasive positive pressure ventilation for stable chronic obstructive pulmonary disease. *Cochrane Database Syst Rev*. 2013;2013(6):Cd002878.
19. Elliott MW. Domiciliary non-invasive ventilation in stable COPD? *Thorax*. 2009;64(7):553-556.
20. Schwarz SB, Magnet FS, Windisch W. Why High-Intensity NPPV is Favourable to Low-Intensity NPPV: Clinical and Physiological Reasons. *Copd*. 2017;14(4):389-395.
21. Dreher M, Storre JH, Schmoor C, Windisch W. High-intensity versus low-intensity non-invasive ventilation in patients with stable hypercapnic COPD: a randomised crossover trial. *Thorax*. 2010;65(4):303-308.
22. Díaz O, Bégin P, Torrealba B, Jover E, Lisboa C. Effects of noninvasive ventilation on lung hyperinflation in stable hypercapnic COPD. *The European respiratory journal*. 2002;20(6):1490-1498.
23. Murphy PB, Brignall K, Moxham J, Polkey MI, Davidson AC, Hart N. High pressure versus high intensity noninvasive ventilation in stable hypercapnic chronic obstructive pulmonary disease: a randomized crossover trial. *Int J Chron Obstruct Pulmon Dis*. 2012;7:811-818.
24. Duiverman ML, Huberts AS, van Eykern LA, Bladder G, Wijkstra PJ. Respiratory muscle activity and patient-ventilator asynchrony during different settings of noninvasive ventilation in stable hypercapnic COPD: does high inspiratory pressure lead to respiratory muscle unloading? *Int J Chron Obstruct Pulmon Dis*. 2017;12:243-257.
25. Murphy PB, Rehal S, Arbane G, et al. Effect of Home Noninvasive Ventilation With Oxygen Therapy vs Oxygen Therapy Alone on Hospital Readmission or Death After an Acute COPD Exacerbation: A Randomized Clinical Trial. *Jama*. 2017;317(21):2177-2186.
26. Struik FM, Sprooten RT, Kerstjens HA, et al. Nocturnal non-invasive ventilation in COPD patients with prolonged hypercapnia after ventilatory support for acute respiratory failure: a randomised, controlled, parallel-group study. *Thorax*. 2014;69(9):826-834.
27. Cherian M, Adam V, Ross B, Bourbeau J, Kaminska M. Mortality in individuals with COPD on long-term home non-invasive ventilation. *Respiratory medicine*. 2023;218:107378.
28. Galli JA, Krahne JS, James Mamary A, Shenoy K, Zhao H, Criner GJ. Home non-invasive ventilation use following acute hypercapnic respiratory failure in COPD. *Respiratory medicine*. 2014;108(5):722-728.
29. Coughlin S, Liang WE, Parthasarathy S. Retrospective Assessment of Home Ventilation to Reduce Rehospitalization in Chronic Obstructive Pulmonary Disease. *J Clin Sleep Med*. 2015;11(6):663-670.
30. Frazier WD, Murphy R, van Eijndhoven E. Non-invasive ventilation at home improves survival and decreases healthcare utilization in medicare beneficiaries with Chronic Obstructive Pulmonary Disease with chronic respiratory failure. *Respiratory medicine*. 2020;177:106291.
31. Frazier W, van Eijndhoven E, Murphy R, Jena AB. Noninvasive ventilation at home reduces mortality in COPD with CRF. *Am J Manag Care*. 2021;27(9):e308-e315.
32. Ergan B, Oczkowski S, Rochweg B, et al. European Respiratory Society guidelines on long-term home non-invasive ventilation for management of COPD. *The European respiratory journal*. 2019;54(3).
33. Macrea M, Oczkowski S, Rochweg B, et al. Long-Term Noninvasive Ventilation in Chronic Stable Hypercapnic Chronic Obstructive Pulmonary Disease. An Official American Thoracic Society

- Clinical Practice Guideline. *American journal of respiratory and critical care medicine*. 2020;202(4):e74-e87.
34. Kaminska M, Rimmer KP, Nonoyama M, et al. Long-term non-invasive ventilation in patients with chronic obstructive pulmonary disease (COPD): 2021 Canadian Thoracic Society Clinical Practice Guideline update. *Can J Respir Crit Care Sleep Med* 2021;5(3):160-183.
  35. Duiverman ML, Vonk JM, Bladder G, et al. Home initiation of chronic non-invasive ventilation in COPD patients with chronic hypercapnic respiratory failure: a randomised controlled trial. *Thorax*. 2020;75(3):244-252.
  36. Masefield S, Vitacca M, Dreher M, et al. Attitudes and preferences of home mechanical ventilation users from four European countries: an ERS/ELF survey. *ERJ Open Res*. 2017;3(2).
  37. Gonzalez J, Sharshar T, Hart N, Chadda K, Raphaël JC, Lofaso F. Air leaks during mechanical ventilation as a cause of persistent hypercapnia in neuromuscular disorders. *Intensive Care Medicine*. 2003;29(4):596-602.
  38. Genta PR, Kaminska M, Edwards BA, et al. The Importance of Mask Selection on Continuous Positive Airway Pressure Outcomes for Obstructive Sleep Apnea. An Official American Thoracic Society Workshop Report. *Ann Am Thorac Soc*. 2020;17(10):1177-1185.
  39. Kaminska M, Montpetit A, Mathieu A, Jobin V, Morisson F, Mayer P. Higher effective oronasal versus nasal continuous positive airway pressure in obstructive sleep apnea: effect of mandibular stabilization. *Canadian respiratory journal*. 2014;21(4):234-238.
  40. Windisch W, Haenel M, Storre JH, Dreher M. High-intensity non-invasive positive pressure ventilation for stable hypercapnic COPD. *Int J Med Sci*. 2009;6(2):72-76.
  41. Orr JE, Coleman J, Criner GJ, et al. Automatic EPAP intelligent volume-assured pressure support is effective in patients with chronic respiratory failure: A randomized trial. *Respirology (Carlton, Vic.)*. 2019;24(12):1204-1211.
  42. Crisafulli E, Manni G, Kidonias M, Trianni L, Clini EM. Subjective sleep quality during average volume assured pressure support (AVAPS) ventilation in patients with hypercapnic COPD: a physiological pilot study. *Lung*. 2009;187(5):299-305.
  43. Magdy DM, Metwally A. Auto-titrating versus fixed-EPAP intelligent volume-assured pressure support (iVAPS) ventilation in patients with COPD and hypercapnic respiratory failure. *Adv Respir Med*. 2021;89(3):277-283.
  44. Gursel G, Zerman A, Basarik B, Gonderen K, Aydogdu M, Memmedova S. Noninvasive auto-titrating ventilation (AVAPS-AE) versus average volume-assured pressure support (AVAPS) ventilation in hypercapnic respiratory failure patients. *Intern Emerg Med*. 2018;13(3):359-365.
  45. Czerwaty K, Dżaman K, Sobczyk KM, Sikorska KI. The Overlap Syndrome of Obstructive Sleep Apnea and Chronic Obstructive Pulmonary Disease: A Systematic Review. *Biomedicines*. 2022;11(1).
  46. Zannin E, Milesi I, Porta R, et al. Effect of nocturnal EPAP titration to abolish tidal expiratory flow limitation in COPD patients with chronic hypercapnia: a randomized, cross-over pilot study. *Respir Res*. 2020;21(1):301.
  47. Windisch W, Vogel M, Sorichter S, et al. Normocapnia during nIPPV in chronic hypercapnic COPD reduces subsequent spontaneous PaCO<sub>2</sub>. *Respiratory medicine*. 2002;96(8):572-579.
  48. Aboussouan LS, Mireles-Cabodevila E. Sleep-Disordered Breathing in Neuromuscular Disease: Diagnostic and Therapeutic Challenges. *Chest*. 2017;152(4):880-892.
  49. Gentile MA. Cycling of the mechanical ventilator breath. *Respir Care*. 2011;56(1):52-60.
  50. Murata S, Yokoyama K, Sakamoto Y, et al. Effects of inspiratory rise time on triggering work load during pressure-support ventilation: a lung model study. *Respir Care*. 2010;55(7):878-884.
  51. Gonzalez-Bermejo J, Rabec C, Perrin C, et al. Treatment failure due to use of a 'ramp' option on pressure support home ventilators. *Thorax*. 2016;71(8):764-766.

- 480 52. Adler D, Perrig S, Takahashi H, et al. Polysomnography in stable COPD under non-invasive  
481 ventilation to reduce patient-ventilator asynchrony and morning breathlessness. *Sleep Breath.*  
482 2012;16(4):1081-1090.
- 483 53. Arnal J-M, Oranger M, Gonzalez-Bermejo J. Monitoring Systems in Home Ventilation. *Journal of*  
484 *Clinical Medicine.* 2023;12(6):2163.
- 485 54. Georges M, Rabec C, Monin E, et al. Monitoring of noninvasive ventilation: comparative analysis  
486 of different strategies. *Respiratory Research.* 2020;21(1):324.
- 487 55. Choi P, Adam V, Zielinski D. Noninvasive Ventilation Downloads and Monitoring. *Sleep Med Clin.*  
488 2020;15(4):569-579.
- 489 56. Toussaint M, Chatwin M, Gonçalves MR, et al. Mouthpiece ventilation in neuromuscular  
490 disorders: Narrative review of technical issues important for clinical success. *Respiratory*  
491 *medicine.* 2021;180:106373.
- 492 57. Majorski DS, Magnet FS, Thilemann S, Schmoor C, Windisch W, Schwarz SB. Portable NIV for  
493 patients with moderate to severe COPD: two randomized crossover trials. *Respir Res.*  
494 2021;22(1):123.

495